

***2nd CRISPR-Based Therapy
Analytical Development Summit 2024***

Analyzing Guide RNA (gRNA) Impurities


Eric Anderson, PhD

18-Sep-2024




We use genome engineering to make healthy cells invisible to drugs


Problem



Few unique cancer antigens, so drugs kill both cancer and healthy cells through **on-target toxicity**



**Vor Paradigm:
Engineered HSCs (eHSCs)**




Cancer cell

Remove target expression on healthy cells so that killing is **cancer-specific**






Trem-cel uses CRISPR/Cas9 editing to delete CD33 in hematopoietic cell transplants during treatment of Acute Myeloid Leukemia






CRISPR/Cas9 effectively and safely edits loss of CD33 from CD34+ cells




Trem-cel has no significant off-target effects






Trem-cel is manufactured at clinical-scale with robust CD33 editing and demonstrates 100% engraftment in 18 patients

Trem-cel clinical-scale manufacturing (30 batches)





Trem-cel demonstrates 100% engraftment



High CD33 editing efficiency (median 89%, range 71-94%)




sgRNA is critical for RNP formation and on-target editing




Identity

```
mGsmA smGsmUCCGA
GCAGAAGAAGAAGUU
UUAGAGCUAGAAAUA
GCAAGUUAAAUAAG
GCUAGUCCGUUAUCA
ACUUGAAAAGUGGC
ACCGAGUCGGUGCmU
smU smUsU
```

Purity




Potency





sgRNA is critical for on-target editing specificity



Safety

	spacer sequence	PAM
sgRNA	GAGUCCGAGCAGAAGAAGAA	NGG
DNA	CTGAGTCCGAGCAGAAGAAGAA	GGGCT
	On-target with perfect matches	

sgRNA	GAGUCCGAGCAGAAGAAGAA	NGG
DNA	TTGA A TCCGAG A AGAAGAAGAA	GGGAC
	Off-target with 2 mismatches	


sgRNA	GA A UCCGAGCAGAAGAAGAA	NGG
DNA	TTGA A TCCGAG A AGAAGAAGAA	GGGAC
	Off-target with 1 mismatch	

sgRNA impurity
G>A substitution

Identity

```
mGsmA smGsmUCCGA
GCAGAAGAAGAAGUU
UUAGAGCUAGAAAUA
GCAAGUUAAAUAAG
GCUAGUCCGUUAUCA
ACUUGAAAAGUGGC
ACCGAGUCGGUGCmU
smU smUsU
```

Purity





What is the FDA guidance for sgRNA as a critical component?



Question-and-Answer



What are FDA's
recommendations for gRNA
purity analysis?

In general, concern is with safety:

Optimize the GE components to reduce the potential for off-target genome modification, to the extent possible.

sgRNA Purity Recommendations

- $\geq 80\%$ purity (HPLC)
- Identify impurities $\geq 1\%$

If $< 80\%$ purity

- Justify purity
- Conduct risk assessment for off-target


sgRNA Sequencing (NGS) Comments

- Not just for identity confirmation
- Orthogonal purity assay to LC/MS




Synthesizing sgRNA at large scale with high purity is challenging

sgRNA synthesis cycle



Coupling efficiency impact on sgRNA full-length product yield




- sgRNA synthesis proceeds from 3' to 5' direction
- 100-mer gRNA \approx 400+ sequential chemical reactions

- Theoretical % FLP = $(\text{efficiency})^{(n-1)}$
- Greater oligo length, less full-length product





Examples of product-related sgRNA impurities







Using NGS and LC/MS as complementary methods for a comprehensive picture of sgRNA purity, identity and impurities

Method	Strengths	Weakness
<p data-bbox="107 546 214 586">NGS</p> 	<ul data-bbox="614 518 1327 732" style="list-style-type: none">• Full length sequence coverage• Impurities sequence / position• High sensitivity (low LLOD)• High specificity	<ul data-bbox="1406 518 2390 732" style="list-style-type: none">• Indirect method• Biases from library prep / sequencing errors• Modifications not captured / detected• Modifications interfere with library prep
<p data-bbox="71 932 224 972">LC-MS</p> 	<ul data-bbox="614 893 1217 1160" style="list-style-type: none">• Direct RNA detection• Intact mass of sgRNA• Impurities identification / quantification• Detection of modifications	<ul data-bbox="1406 893 2351 1051" style="list-style-type: none">• Intact mass alone not sufficient for identity• sgRNA sequence confirmation• Impurity sequence / position



Using NGS and LC/MS as complementary methods for a comprehensive picture of sgRNA purity, identity and impurities

Purity





Purity as integrated peak area

Challenges:

- Peak tailing / overlapping
- Single-nucleotide resolution

Identity





LC-MS: Intact mass



NGS: Full-length sequencing

LC-MS/MS: Full-length sequencing (challenging)

NGS



Impurities




LC/MS: identification, quantification, chemical modifications / impurities

NGS: position of sequence impurities




HPLC is used for routine sgRNA purity assessment, but single-nucleotide resolution is challenging



100nt sgRNA with impurity spike ins

Green = n+1
 Black = n-1
 Pink = n-2
 Brown = n-5
 Blue = FLP + Imp.





Challenges:

- Peak tailing leads to peak overlapping
- Resolving and collecting impurities
- Applying consistent peak integration



LC-MS/MS sequencing method provides complete coverage and sequence confirmation of modified sgRNA



Full Sequence Confirmation

- LC/MS detected both full length and partial 5' and 3' cleavage products
- MS/MS sequenced chemically modified 5' and 3' ends
- Complete coverage provided by overlapping sequencing of 5' and 3' cleavage products
- Assay specific for test sgRNA – non-target sgRNA or n-1 fail sequence identity
- Submitted in VOR33/trem-cel IND


Confidential

MS-Based Sequencing Method Presented by Vor Bio and Axolabs at 2021 TIDES






NGS confirms sequence identity of sgRNA



Alignment	Mismatch+gap	Replicate 1			Replicate 2			% reads aligned to reference sequence
		Test Sample 1	Test Sample 2	Test Sample 3	Test Sample 1	Test Sample 2	Test Sample 3	
XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAAAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCAGTCGGTGCT XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAAAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCAGTCGGTGCT	0	65.8	61.7	71.8	63.7	59.3	69.1	




Full Sequence Confirmation



- same sgRNA sequence in 3 different batches, 2 replicates
- 3' terminal UUUU trimmed for alignment
- most frequent sequence matches expected reference sequence
- >1M reads support full-length alignment
- sgRNA purity influences % reads aligned to reference sequence



Over 80% of NGS reads account for the full length and n-1 sgRNA sequences while all other sequences individually are <1%



Percentage of reads mapped to sgRNA full length and impurity sequences




Sequence	Alignment	Mismatch+gap	Test Sample 3 (read counts)	Test Sample 3 (% reads)
FLP	XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT	0	1.77E+06	71.75
n-1	XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT ~XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT -XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT	1	3.05E+05	12.38
n-2	XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT ~XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT --XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT	2	4253	0.17
internal deletion	XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTG-AAAAAGTGGCACCGAGTCGGTGCT	1	10334	0.42
n+1	-XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT ~XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT	1	9873	0.40
internal deletion	XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTCGGTGCT XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTC- XXXXXXXXXXXXXXXXXXXXGTTTGTAGAGCTAGAAATAGCAAGTTAAAATAAGGCTAGTCCGTTATCAACTTGAAAAAGTGGCACCGAGTC-GTGCT	1	1578	0.06




NGS reveals majority of insertions and deletions in sgRNA are <1% and tend to occur at the 5' end

Insertions and deletions in the sgRNA sequence




Cumulative frequency of sgRNA deletions by position






sgRNA sequence impurities identified by NGS are likely a combination of actual impurities and errors during library prep and sequencing






LC/MS fragment analysis allows for analysis of sgRNA chemical modifications and side reaction-derived impurities



5' cleavage fragment (27nt) with impurity spike ins





Assay design / optimization:

- Focused on LC/MS analysis of sgRNA spacer by generating short 5' cleavage fragment (27nt)
- RNaseH / ASO design, conditions, activity, specificity
- IP-RP-HPLC method
 - no RNaseH control / blank subtraction
 - 10% truncated impurity spike ins
 - resolves n-1, n-2, n-5 truncations
 - note truncated series peaks (e.g., n-1) <<10%
 - does not fully resolve n+1



LC/MS fragment analysis confirms high purity and identifies impurities within the spacer region of sgRNA

5' cleavage fragment (27nt) with mass assignments



5' cleavage fragment (27nt) analysis:


- precise cutting of 27nt, 5' fragment
- main peak is 27nt (88%)
- Main peak observed MW matches theoretical MW of FLP

Impurities:

- Late-eluting peaks include n+1 and methylamine adducts
- truncation series n-1, n-2, etc.
- n-1 is the terminal 2'OMe base (<1%)
- abasic sites observed
- phosphorothioate impurities (PS->PO)
- degradation (e.g., 3'(n-1) +2'-3'-cyclic-P)



Overall sgRNA impurity strategy – how comprehensive is “comprehensive characterization”?



sgRNA Impurity Considerations

- What is the sgRNA purity?
- Can impurities be detected, identified and quantified?
- What type of impurities and where are they in the sgRNA sequence?
- What is the potential off-target risk and is this a novel off-target?
- What are the impacts on RNP formation, on-target editing and manufacturing process performance?
- Is isolation or synthesis of impurities needed for further characterization and functional testing?



Summary

- sgRNA is a critical component required for genome engineering
- GMP manufacturing of sgRNA at scale is challenging, resulting in product-related impurities that may require identification
- sgRNA purity has the potential to impact both on and off-target editing activity
- sgRNA impurity identification is a regulatory requirement and can improve manufacturing processes
- HPLC, NGS and LC-MS/MS methods are complementary in assessing sgRNA identity, purity and impurity identification



Acknowledgements

Vor

Tirtha Chakraborty

John Lydeard

Gary Ge

Jianxin Hu

Lakshmi Karthik

Antonino Montalbano

Huan Qiu

Ruijia Wang

Sawyer Letourneau

Axolabs

Madeleine Mohrand

Jasmin Brugger

Julia Paulus

Stephan Seiffert



www.vorbio.com

Confidential